



LivaNova

Advanced Circulatory Support

LifeSPARC Circuit Emergency Guide

24/7 Clinical Support 800 373 1607

Customer Service 800 373 7421

Office 620 Alpha Dr, Suite 2
Pittsburgh, PA 15238

Visit LivaNova.com

IM-7100219-ACS

Controller Critical Failure**Possible
Presentation**

- Continuous audible controller alarm
- Critical failure icon on main display

**Immediate
Action**

- **Do not turn off controller**
- **DO NOT DISCONNECT OR STOP THE PUMP**
unless determination has been made to discontinue support
- Confirm that the pump is running and that the pump speed is being maintained as displayed on the secondary display LED lights
- Controller replacement will be necessary to re-establish full functionality
- Contact LivaNova for additional guidance and support
 - Local representative
 - **24/7 CLINICAL SUPPORT:**
800-373-1607

Next Steps

- Notify physician in charge and appropriate team members
- Perform LifeSPARC controller exchange
- *Refer to LifeSPARC Operations Manual for detailed instructions*

Pump Failure**Possible
Presentation**

- Rule out patient factors, ensure controller is on and functioning appropriately
- Potential controller alarms
 - “Pump stopped”, “Low Flow”, “Pump current high”
- Patient deterioration
- Low flow or no flow through the circuit

Immediate Action

- Clamp both inflow & outflow limbs of circuit
- Support patient
 - Hemodynamic & ventilator support as needed
- Notify physician in charge and appropriate team members

Next Steps

- Prepare for pump or circuit exchange
- Contact LivaNova for additional guidance and support
 - Local representative
 - **24/7 CLINICAL SUPPORT:**
800-373-1607

Unintentional Decannulation**Possible Presentation**

- Rapid decline in flows and patient deterioration
- Massive hemorrhage or rapid swelling & discoloration at cannulation site
- Cannula visibly out of insertion site

Immediate Action

- “3C rule”
 - **CLAMP** both inflow & outflow limbs of circuit, stop flow & pump
 - **CALL** for help
 - **COMPRESS** firmly proximal to cannula insertion site

Next Steps

- Support patient
 - Hemodynamic & ventilator support as needed
 - Massive transfusion
- Notify physician in charge and appropriate team members

Mitigations

- Verify appropriate placement and securement of all cannulas regularly and prior to mobilization
- Monitor and record insertion depth of all cannulas
- Verify patient bed and circuit are secure and in appropriate position

Erratic Flows / Suction Event**Possible Presentation**

- Sudden, large fluctuations in flows despite no change in RPMs
- (+/-) chattering present in venous drainage line
- Patient deterioration

Immediate Action

- Attempt to stabilize flows – decrease RPMs until stable flow achieved
- Complete circuit check: venous drainage cannula position, kinks, bleeding, thrombus
- Notify physician in charge and appropriate team members
- Consider volume challenge if hypovolemia suspected

Next Steps

- Identify and treat underlying cause
 - Potential patient causes: hypovolemia, hemorrhage, tamponade, excessive intraabdominal or intrathoracic pressure
 - Potential circuit causes: cannula malposition, circuit kinking, unnecessarily high RPMs, or clot in circuit

Mitigations

- Ensure adequate patient volume status
- Rapid diagnosis and treatment of hemorrhage, tamponade, excessive intraabdominal, intrathoracic pressure
- Regular monitoring of circuit and verification of cannula placement

Air Entrainment**Possible Presentation**

- Air trapped in top of oxygenator
- Visible air in tubing or pump head
- No flows due to air-locked pump

Immediate Action

- Clamp both inflow & outflow limbs of circuit
- Place patient in Trendelenburg position if air delivered to patient
- Support patient
 - Hemodynamic & ventilator support as needed
- Notify physician in charge and appropriate team members

Next Steps

- De-air the circuit or complete circuit exchange

Mitigations

- Minimize number of circuit access sites
- Ensure all connections proximal to pump are occlusive and secure
- Verify appropriate placement and securement of all cannulas regularly
- Ensure proper care when accessing central lines

Cardiac Arrest

Possible Presentation	<ul style="list-style-type: none">• Cardiac arrest or arrhythmia
Immediate Action	<ul style="list-style-type: none">• VV ECMO:<ul style="list-style-type: none">• Initiate CPR• VA ECMO:<ul style="list-style-type: none">• Maximize ECMO pump flow (may require volume)• Initiate CPR if unable to maintain adequate flow
Next Steps	<ul style="list-style-type: none">• Evaluate for and treat underlying cause of cardiac arrest
Mitigations	<ul style="list-style-type: none">• Monitor for and treat electrolyte imbalances, acid base abnormalities, hypoxia, bleeding/hypovolemia, cardiac tamponade, tension pneumothorax, etc.



Manufactured by CardiacAssist, Inc
620 Alpha Drive, Pittsburgh, PA 15238, USA
412-963-7770

24/7 Clinical Support 800 373 1607

Customer Service 800 373 7421

Office 620 Alpha Dr, Suite 2
Pittsburgh, PA 15238

Visit LivaNova.com

IM-7100219-ACS